





# Container Packaging and Steam Sterilization

BY TAYLOR HAMSCHER, ASSOCIATE PRODUCT MANAGER, HEALTHCARE SOLUTIONS—AESCULAP INC.

Certified Instrument Specialist (CIS) lessons provide members with ongoing education in the complex and ever-changing area of surgical instrument care and handling. These lessons are designed for CIS technicians, but can be of value to any CRCST technician who works with surgical instrumentation.

## Earn Continuing Education Credits

**Online:** Visit [www.myhspa.org](http://www.myhspa.org) for online grading.

**By mail:** For written grading of individual lessons, send completed quiz and \$15 to: Purdue Online – Central Service Young Hall, Room 405 155 S. Grant St. West Lafayette, IN 47907

**Subscription series:** Purdue Extended Campus offers an annual mail-in or online self-study lesson subscription for \$75 (six specific lessons worth two credits each toward CIS recertification of six credits). Call 765.494.0959 for details.

**Scoring:** Each quiz with a passing score is worth two continuing education (CE) credits toward your CRCST recertification (12 credits).

**More information:** HSPA provides online grading service for any of the Lesson Plan varieties. **Purdue University provides grading services ONLY for CRCST and CIS lessons.** Direct questions about online grading to HSPA at 312.440.0078. Questions about written grading are answered by Purdue University at 765.494.0959.

## LEARNING OBJECTIVES

1. Review basic packaging requirements
2. Identify how rigid containers react during the steam sterilization process
3. Explore characteristics of rigid containers

**M**aintaining sterility of items used in the Operating Room (OR) is essential for patient safety. Multiple types of packaging are available to help maintain sterility. This lesson addresses one type—rigid containers—and their sterilization process.

### Objective 1: Review basic packaging requirements

There are multiple sterilization modalities available for Sterile Processing (SP), including steam, low-temperature methods and, occasionally, ethylene oxide (EO). Regardless of the sterilization method, the type of packaging material used must comply with key performance features. These include allowing penetration and evacuation of the sterilant; maintaining sterility of sterilized contents until they are opened; allowing for aseptic presentation; and having the ability to label the packaging.

Three primary packaging methods are offered today: paper/plastic peel pouches, flat wrappers, and rigid sterilization containers. When evaluating selecting packaging methods, it is important to ensure that the packaging/container:

- Offers protection during transport, storage and handling
- Is free of toxic ingredients or non-fast dyes
- Has minimal lint
- Can close securely (with a tamper-evident seal)
- Resists tears and punctures
- Is correctly sized to ensure even distribution of items to be sterilized
- Offers ease of use when preparing and opening the packaging
- Includes the manufacturer's instruction for use (IFU)
- Is cost effective

Paper/plastic peel pouches are designed for small, often individual, and



lightweight instruments. The sterilant enters the pouch, is removed on the paper side and is placed in the sterilizer to promote proper sterilant access. Labeling should be done on the plastic side of the pouch. Flat wrappers come in a wide range of sizes to cover the devices being sterilized. It is important to use properly sized packaging; too-large wrappers will have excess material that could cause steam to accumulate and pool, causing a wet pack. Proper wrapping technique ensures that the first fold completely covers the package and that the wrapper remains snug throughout the wrapping process. The packaging should be inspected each time it is handled to ensure no holes or tears are present.

Rigid sterilization containers are most commonly used when sterilizing multiple items (often, instrument sets). Containers come in a variety of sizes and feature filters or valves for sterilant penetration. Filters are often made of a cellulose or polypropylene material, both of which allow for penetration and evacuation of the sterilant and provide a long-term sterile barrier during storage. For sterilization containers, it is important to ensure the combined weight of the container and contents does not exceed 25 pounds (per Association for the Advancement of Medical Instrumentation standards and Association of periOperative Registered Nurses guidelines). They must also be easy to carry and handle; have prepping features such as loading, closing and locking; have stacking capabilities for transport and shelf storage (as well as stacking capability for internal baskets); feature effective filter retention mechanisms; come in a variety of sizes; and be compatible with various sterilization modalities.

For all types of packaging material, it is critical to understand the shelf life for sterilized items. Two considerations

for shelf life are expiration date and event-related practice. Regardless of the shelf-life method a facility employs, it is crucial to always follow industry and manufacturer guidelines and IFU. Standards ANSI/AAMI ST77 and ST79 provide recommendations for rigid containers and steam sterility and assurance, respectively. Be sure to check with each packaging manufacturer to determine if there is an expiration date assigned to the packaging material.

Additionally, when ensuring effective packaging to maintain sterility, aseptic presentation must also be understood and performed properly. Aseptic presentation involves transferring sterile contents from their sterile barrier packaging to a sterile field using procedures that minimize the risk of microbial contamination. *Note: This is a critical step in the process to keep patients safe. Remember to always follow the packaging manufacturer's IFU, examine and inspect trays, baskets, liners, containers, etc., ensure labeling and indicators are easy to see and read, and check for packaging integrity after sterilization and before issuing the device.*

### Objective 2: Learn how rigid containers react during the steam sterilization process

Prior to assembly and sterilization it is important to properly clean, rinse and dry the container following the manufacturer's IFU. Rigid containers must be cleaned after each use, even when they do not visibly appear soiled. During steam sterilization, items inside the containers are heated to the manufacturer's recommended appropriate temperature for the correct period amount of time. Rigid containers offer several advantages, including greater protection of contents, durability, cost-effectiveness, and increased staff efficiency and throughput.

Once loaded, sealed, locked and labeled, the rigid container is ready to be placed in the sterilizer. A pre-vacuum sterilizer goes through three phases in a cycle: conditioning, exposure and exhaust. During the conditioning phase, the dynamic air removal process removes air, and the chamber is heated. The chamber continues to heat until the set sterilization temperature is reached. During this time, steam penetrates the container filter, surrounding and warming the contents. In the exposure phase, the chamber holds the set temperature for the required amount of time to sterilize the chamber contents. In the final step, the exhaust phase, air is removed from the chamber and pressure is gradually reduced. Content drying also occurs during this phase. Upon completion of the three phases, the chamber can be opened, and the cooling process can begin. *Note: Do not touch the containers until the contents have cooled to room temperature.* The filter acts as a part of the sterile barrier, preventing contamination of the tray contents.



Figure 1: Placing filter in the container lid





### Objective 3: Explore characteristics of rigid containers

The Association for the Advancement of Medical Instrumentation (AAMI) notes that rigid containers are “valuable tools in sterility maintenance and to prevent damage to surgical instruments and devices.” Rigid sterilization container systems are durable and comprised of a rigid material such as anodized aluminum, stainless steel, or proprietary polymers.

These containers are validated by the U.S. Food and Drug Administration as Class II medical devices, which indicates the container has been designed, tested and validated to conform to industry standards in terms of specified sterilization modality, transport/shelf life, aseptic presentation, and labeling. The durable material used to make the containers ensures that the containers can endure multiple sterilization cycles, handling, storage, etc.



Figure 2: Outer container with inner trays

When properly maintained, many hospitals and ambulatory surgery centers can last 10 years or more, maximizing the return on investment. Proper care and handling that follows the manufacturer’s IFU is paramount, along with inspection before and after each use.

### Conclusion

Approved rigid containers have been FDA approved through the 510(k) process and have been on the market for more than 40 years, representing a lengthy history of use and efficacy in U.S. healthcare facilities. When properly handled and maintained, rigid containers can prevent damage to surgical devices, maintain the sterility of their contents until the processed devices are used in a procedure, and provide many years of use. **P**

### RESOURCES

1. Klacik S. “Proper packaging ensures integrity of sterilization. *OR Manager*. Feb. 19, 2016.
2. Seavey R. “Packaging for Sterilization: It’s Not Just a Box, a Bag or Gift Wrapping.” 3M Health Care lesson. *Managing Infection Control*. Jan. 2008.
3. Meggs J. “Thinking of Buying ... Rigid Sterilization Containers.” *Outpatient Surgery*. June 2, 2016.
4. STERIS. “Steam Sterilization for Medical Equipment.” Aug. 16, 2017. [www.steris.com/healthcare/knowledge-center/sterile-processing/steam-sterilization-for-medical-equipment](http://www.steris.com/healthcare/knowledge-center/sterile-processing/steam-sterilization-for-medical-equipment)
5. Aesculap USA. “SterilContainer™ System Performance Overview.” Nov. 1, 2016.
6. Censis. “Rigid Containers for Sterilization.” <https://censis.com/blog/rigid-containers-for-sterilization/>
7. Ortiz R. “The Benefits of Rigid Sterilization Containers.” *Outpatient Surgery*. Aug. 2, 2022.
8. Practice Greenhealth. “Rigid Container Systems.” <https://practicegreenhealth.org/topics/greening-operating-room/rigid-sterilization-containers>



# CIS Self-Study Lesson Plan Quiz

## Container Packaging and Steam Sterilization

Lesson No. CIS 297 (Instrument Continuing Education – ICE) · Lesson expires May 2026

1. Which of the following is an approved, effective packaging system?
  - a. Peel pouches
  - b. Flat wrapper
  - c. Rigid sterilization containers
  - d. All the above
2. Which is a key performance feature for effective packaging?
  - a. It allows for penetration and evacuation of the sterilant
  - b. It maintains sterility of the contents until opened
  - c. It allows for aseptic presentation
  - d. All the above
3. Paper/plastic peel pouches:
  - a. Are designed for small groups of instruments or single, heavier instruments
  - b. Allow the sterilant to enter and be removed on the plastic side
  - c. Require that labeling be done on the plastic side
  - d. Are used to contain like devices inside an instrument set
4. Flat wrappers:
  - a. Come in a wide range of sizes
  - b. Can be stacked in sterile storage
  - c. Can be reused as long as the tray is unused
  - d. Are rarely used in healthcare
5. Rigid sterilization containers:
  - a. Are most commonly used for small sets and delicate devices
  - b. Use filters or valves
  - c. Have no weight restrictions
  - d. Are the least cost-effective type of packaging
6. Sterilization containers should not exceed 20 pounds, including the container and its contents.
  - a. True
  - b. False
7. Which Association for the Advancement of Medical Instrumentation (AAMI) standard provides recommendations for rigid containers?
  - a. ST59
  - b. ST79
  - c. ST32
  - d. ST11
8. Aseptic presentation:
  - a. Is the process of transferring sterile contents from sterile barrier packaging
  - b. Involves effectively sealing the package for sterilization
  - c. Involves the entire cleaning and sterilization process
  - d. Is only used in the Operating Room
9. Which type of packaging material can be sterilized in a steam sterilizer?
  - a. Rigid container
  - b. Peel pouch
  - c. Flat wrap
  - d. All the above
10. What is the first phase of steam sterilization?
  - a. Conditioning
  - b. Cooling
  - c. Pre-vacuum exhaust
  - d. None of the above
11. A primary advantage of using rigid containers is:
  - a. They can be stacked as high as necessary due to their durability
  - b. No instructions for use are required
  - c. Greater protection of contents
  - d. They can be used indefinitely
12. The U.S. Food and Drug Administration (FDA) labels rigid containers as:
  - a. Class I medical devices
  - b. Class II medical devices
  - c. Class III medical devices
  - d. The FDA does not classify rigid containers
13. When selecting a packaging system, which of the following should be considered in the decision-making process?
  - a. Container shape
  - b. ANSI/AAMI ST91 recommendations
  - c. Container weight
  - d. FDA guidelines
14. The shelf life of containerized items is based on which of the following?
  - a. Number of items in the container
  - b. Manufacturer's instructions for use
  - c. Type of sterilization method used
  - d. All the above
15. After sterilization, containers should cool:
  - a. For at least one hour
  - b. For two hours
  - c. Until the contents reach room temperature
  - d. Until the contents feel cool enough to handle

### Request for paper/pencil scoring

(Please print or type information)

- I have enclosed the scoring fee of \$15. (Please make checks payable to Purdue University. We regret that no refunds can be given.)
- Check here if you have a change of address.
- Check here if you wish to have your results emailed to you.

### Detach quiz, fold and return to:

Purdue Online – Central Service  
 Young Hall Room 405  
 155 S. Grant St.  
 West Lafayette, IN 47907

Name \_\_\_\_\_

Mailing Address (Be sure to include apartment number or post office box.) \_\_\_\_\_

City \_\_\_\_\_

State/Province \_\_\_\_\_

Zip Code/Postal Code \_\_\_\_\_

Daytime Telephone \_\_\_\_\_

HSPA Membership Number \_\_\_\_\_

Email Address \_\_\_\_\_

*If your name has changed in the last 12 months, please provide your former name. Purdue University is an equal access/equal opportunity institution.*