

Pass



Fail





Chemical Indicator Classifications: What Technicians Need to Know

CRAIG WALLACE, PRESIDENT—WALLACE STERILIZATION CONSULTING

This series of self-study lessons on SP topics was developed by the HSPA. Certified Registered Central Service Technician (CRCST) lessons provide members with Purdue University's Extended Campus and HSPA both offer grading opportunities.

Earn Continuing Education Credits

Online: Visit www.myhspa.org for online grading.

By mail: For written grading of individual lessons, send completed quiz and \$15 to: Purdue Online – Central Service Young Hall, Room 405 155 S. Grant St. West Lafayette IN, 47907

Scoring: Each quiz with a passing score is worth 2 contact hours toward your CRCST recertification (12 hours).

Subscription series: Purdue Extended Campus offers an annual mail-in or online self-study lesson subscription for \$75 (six specific lessons worth 2 points each toward CRCST recertification of 12 hours). Call 765.494.0959 for details.

More information: HSPA provides online grading service for any of the Lesson Plan varieties. **Purdue University provides grading services ONLY for CRCST and CIS lessons.** Direct any questions about online grading to HSPA at 312.440.0078. Questions about written grading are answered by Purdue University at 765.494.0959.

LEARNING OBJECTIVES

1. Review chemical indicator designs and applications
2. Discuss the International Organization for Standardization's chemical indicator classifications
3. Examine the U.S. Food and Drug Administration's approach to chemical indicators and compare it to the International Organization for Standardization's approach

Chemical indicators (CIs) are an important part of a sterilization quality control testing program. CIs are used in conjunction with physical monitors and biological indicators (BIs) to provide critical information on the quality of a sterilization process. This information is necessary to ensure that all instruments are safe and ready for patient use.

Objective 1: Review chemical indicator designs and applications

CIs are defined as a “test system that reveals change in one or more pre-specified process variables based on a chemical or physical change resulting from exposure to a process.”¹ CIs respond to certain specified critical variables in the sterilization process, and their response is chemical or physical as opposed to biological (such as a BI response). There are two main types of CI technology. The first is a reactive

ink chemistry that changes color after exposure to the process (e.g., cream-colored ink darkens after exposure). The second is a chemical tablet enclosed within the indicator that melts under certain process conditions, and the liquid chemical then travels down a window to reach a pass/fail line. Some designs incorporate both concepts: a sterilant penetrates a restricted access and causes a color change as the sterilant travels down the window.

CIs have three primary applications. The first are external or exposure indicators. These color-change indicators are visible from the outside of a package or container and are intended to provide fast visual evidence that the package has been exposed to the sterilization process. The second applications are special tests such as Bowie-Dick test packs. The third applications are pack or internal indicators. These CIs are placed inside packages or containers and



provide information on the physical conditions inside the package, where the instruments are located. The CI performance requirements are different for each application.

CIs are segmented into different classifications to better describe their specific applications and performance requirements. The three primary CI applications (external indicators, special tests and internal indicators) are quite different and require different CI performance characteristics. CIs are classified by two different organizations: the International Organization for Standardization (ISO) and the U.S. Food and Drug Administration (FDA). These organizations have different charters and objectives and classify CIs differently. It is useful to understand the two organizations' CI classifications, so it is easier to understand CI labeling and instructions for use (IFU).

Objective 2: Discuss the International Organization for Standardization's chemical indicator classifications

ISO is dedicated to standardization to promote global trade and equitable financial growth. ISO connects the national standards bodies of 165 countries, and ISO's standardization work is carried out by approximately 265 technical committees that cover subject areas ranging from dentistry to freight containers to photography. ISO Technical Committee 198 (TC198) is focused on the sterilization of healthcare products. Within this Technical Committee are 13 working groups comprised of international experts who develop standards for specific subject areas (such as moist heat sterilization and sterilization packaging). ISO TC198 Working Group 6 is responsible for all international standards related to CIs.² The

| Type | Name | Description |
|------|---|--|
| 1 | Process indicators | Provide visual evidence of exposure to the process |
| 2 | Indicators for specific tests | Special tests, such as Bowie-Dick tests |
| 3 | Single critical process variable indicators | Respond to one critical process variable |
| 4 | Multicritical process variable indicators | Respond to two or more critical process variables |
| 5 | Integrating indicators | Respond to all critical process variables. Stated values must be related to theoretical biological indicator. Multiple stated values are required across the entire sterilization temperature range. |
| 6 | Emulating indicators | Respond to all critical process variables. Single stated value based on the specified cycle. |

Figure 1: ISO chemical indicator types³

standard that establishes the ISO CI classifications is ISO 11140-1.³

The U.S. national version of this standard is ANSI/AAMI/ISO 11140-1; this version is identical in content to ISO 11140-1 and describes labeling and performance requirements for CIs. It is generally considered a manufacturer's standard, but familiarity with its content is useful for Sterile Processing (SP) professionals to help understand CI labeling and the performance capabilities of each type of CI.

CIs that meet the requirements of the ISO standard are typically labeled as compliant with the standard. Compliant CIs should be used whenever possible, because they provide a level of assurance that the chosen indicator's performance is appropriate and based on international norms. It should be noted, however, that CI manufacturers self-certify compliance to the standard; there is no international regulatory organization that independently verifies compliance. Some manufacturers will choose to hire a third-party organization (called a "notified body") to verify their compliance claims.

ANSI/AAMI/ISO 11140-1 establishes six types of CIs. They are numbered Types 1 through 6, but there is no hierarchical implication to the

numbering system (i.e., six is not better than one; they are simply different categories). The six types of CIs and their descriptions are provided in **Figure 1**.

It should be noted that Type 3 indicators still appear in the standard but are not commonly used because the other internal indicators (Types 4, 5 and 6) provide more information and are readily available. It should also be noted that not all indicator types are described or specified for all sterilization processes.

The ISO standard describes the three main applications for CIs (as previously addressed in this lesson). That is, CIs may be used to visually differentiate processed and unprocessed items, they may be used in special tests such as the Bowie-Dick, or they may be used inside individual load items to evaluate process variables and parameters at the point of placement.³

Each type of CI has performance requirements specified in the standard. Learning some CI terminology can help technicians understand how the performance of each type is specified in the standard. The first important term is "stated value" (SV), which is defined as "value or values of a critical variable at which an indicator is designed to reach its endpoint."¹ Reaching its end point



| Exposure Time | Temperature | Expected Result |
|---------------|-------------|-----------------|
| 2 minutes | 121°C | Fail |
| 10 minutes | 121°C | Pass |
| 0.3 minutes | 134°C | Fail |
| 2 minutes | 134°C | Pass |
| 30 minutes | 140°C | Fail |

Figure 2: Type 1 process indicator requirements (steam)³

| Type | Test Point | Time | Temperature |
|------|------------|------|-------------|
| 4 | 1 | SV | SV |
| 4 | 2 | -25% | -2°C |
| 5 | 1 | SV | SV |
| 5 | 2 | -15% | -1°C |
| 6 | 1 | SV | SV |
| 6 | 2 | -6% | -1°C |

Figure 3: Internal indicator performance requirements (steam)³ In Figure 3, test point 1 is the pass condition, and test point 2 is the fail condition. Fail conditions are created by reducing the critical variable values to reduce the effectiveness of the process and create a process failure.

means showing a “pass” result, so the SV is the sterilization exposure conditions where the indicator must show a pass result. For example, if a steam CI has an SV of 132°C for three minutes, it means the CI must show a pass result when exposed to saturated steam at 132°C for three minutes in a laboratory test sterilizer. The SV(s) for a CI must be listed on the product labeling to be compliant with the standard.

Another important term is “critical process variable,” which identifies the elements of a sterilization process that have direct impact on the lethality (or killing power) of the process. For example, the critical process variables for steam sterilization are exposure time, temperature and steam saturation. Any change in these variables (e.g., reduced exposure time or temperature) will affect the process. An example of a non-critical variable for steam sterilization is

pressure. Pressure change alone will not affect the lethality of the process.

In general, the ISO standard defines the sterilization conditions where the indicator must show a pass result (the CI reaches its end-point color or travels past the pass/fail line) and conditions where the indicator must show a fail result. In this way, the standard defines the level of challenge the indicator must present to the sterilization process. *Note: Type 2-specific test indicators have a separate set of ISO standards that specify their requirement.*

There is a difference in how performance requirements are specified for Type 1 process indicators (external indicators to verify exposure) and Type 3, 4, 5 and 6 indicators (intended for use inside packages). Type 1 indicators have a prescribed set of pass/fail conditions, while the internal indicators require the use of SV. This is because Type 1

indicators are only intended to verify exposure to the process and differentiate processed from unprocessed items, whereas internal indicators are expected to provide information on the quality of the process by evaluating the critical process variables.

Figure 2 provides an example of how process indicator performance is specified for steam sterilization. For comparison, Figure 3 shows how steam internal indicator performance is specified.

Note: Since Type 1 indicators are not intended to evaluate the quality of the process, there is no link to SVs in their requirements.

Internal indicators evaluate the quality of the sterilization process by responding to the critical process variables. Type 3 and 4 indicators are only required to respond to one variable or one or more variables, respectively. Type 5 integrating indicators for steam sterilization are required to respond to all of the critical process variables and must also have three SVs specified across the range of sterilization process temperatures. This means that they must have a SV at 121°C, one at 135°C and at least one at a temperature in between (if three SVs are used, the third is usually at 128°C). These SVs should be at least as challenging as the performance of a theoretical BI across the full sterilization range. Type 6 emulating indicators are also required to respond to all critical process variables, but they only have a single SV that must reflect the specific process in which it is used. In other words, Type 6 indicators are cycle specific.

It should be noted that ISO 11140-1 does not define or provide performance requirements for Type 5 and 6 CIs for vaporized hydrogen peroxide (VH₂O₂) sterilization processes. Type 4 multicritical process variable indicators,



therefore, provide the highest level of information for VH2O2 processes.

Objective 3: Examine the U.S. Food and Drug Administration's approach to chemical indicators and compare it to the International Organization for Standardization's approach

The FDA is a federal organization with broad regulatory responsibilities focused on protecting the health and safety of the public. As part of its responsibility, the FDA regulates medical devices. CIs are considered Class II medical devices and are subject to FDA regulation.

This regulation includes a requirement that a CI manufacturer submit a premarket submission for FDA review before the manufacturer can sell the CI to healthcare facilities. This notification document is called a "510(k)" (named for the clause in the Food and Drug Act that defines this process) and contains testing data and labeling information that will be reviewed by the FDA. If the 510(k) document supports the CI as safe and effective, the FDA will "clear" it and allow for commercial distribution. The FDA uses product codes to categorize medical devices for their internal systems, including 510(k). For example, the product code for BIs is "FRC," and the product code for standard CIs is "JOJ."

The FDA published a guidance document for CI manufacturers that outlines the information and testing data that must be included in the 510(k) submission.⁴ This document lists the three CI classifications used by the FDA. These classifications and their descriptions and applications are shown in **Figure 4**.

The FDA uses only three CI classifications, while ISO uses six. Both

| Classification | Description | Intended use |
|---------------------|---|--|
| Process indicator | A chemical indicator intended for use with individual units (e.g., packs, containers) to demonstrate that the unit has been exposed to the sterilization process and to distinguish between processed and unprocessed units | This process indicator is intended to be used by a healthcare provider with sterilization wraps, containers, cassettes or pouches to distinguish between processed and unprocessed units. |
| Chemical integrator | A chemical indicator designed to react to all critical parameters over a specified range of sterilization cycles | This chemical integrator is intended to be used by a healthcare provider to demonstrate that the parameters over a specified range of sterilization cycles have been met in a specified sterilization wrap, container, cassette or pouch. |
| Air removal test | A chemical indicator to be used in the standard test pack to determine the efficacy of the air removal phase in the steam sterilization process | This air removal indicator is intended to be used by a healthcare provider in a standardized test pack (e.g., Bowie-Dick test pack) to assess the ability of a pre-vacuum sterilizer to remove air and allow steam to penetrate into wrapped goods and porous loads. |

Figure 4: CI descriptions and intended use⁴

| Steam | Temperature | Time to show a pass result |
|---------------|-------------|----------------------------|
| Yes | 121°C | 2 minutes to 10 minutes |
| Yes | 132°C–135°C | 20 seconds to 2 minutes |
| No (dry heat) | 140°C | No pass at 30 minutes |

Figure 5: FDA steam process indicator requirements

| Steam | Exposure Time | Temperature | Expected Result |
|---------------|---------------|-------------|-----------------|
| Yes | 2 minutes | 121°C | Fail |
| Yes | 10 minutes | 121°C | Pass |
| Yes | 0.3 minutes | 134°C | Fail |
| Yes | 2 minutes | 134°C | Pass |
| No (Dry heat) | 30 minutes | 140°C | Fail |

Figure 6: ISO steam process indicator requirements

organizations have a classification for process indicators and air removal (special) tests. The FDA places all internal indicators (ISO Types 3, 4, 5 and 6) into a single classification: CIs. The FDA's test requirements are different than the ISO's. ISO requirements tend to be more

specific as the test conditions for both pass and fail tests are defined. Examples of the FDA requirements for steam process indicators are provided in **Figure 5**.

As shown in some of the figures, the ISO requirements are more detailed than the FDA requirements, with



specific conditions defined for both pass and fail results. The FDA requires pass/fail results in the 510(k) document, but the agency leaves it to the manufacturer to determine the test parameters.

For steam integrators, the FDA definition of responding to all critical variables will only compare to ISO Type 5 and Type 6 indicators. In this case, the ISO requirements again are more prescriptive. The FDA testing for CIs asks for side-by-side testing with a BI, while ISO Type 5 indicators require SVs that are related to a calculated BI response.

Requirements for CIs for VH₂O₂ sterilizers are covered under “New Technology Sterilizers” in the FDA guidance document. The FDA does not provide any specific performance criteria for CIs for VH₂O₂ processes but indicates that the general approach used for steam CIs should be applied. The FDA has also accommodated innovative technologies in this area and recently added product code “QKM” for chemical vapor sterilization multivariable CIs. These indicators monitor one or more critical process variables, such as H₂O₂ concentration, temperature and exposure time. This new category aligns closely with the ISO Type 4 definition.

Conclusion

CI classifications are useful as they define the capabilities and appropriate use of the wide range of CIs available. ISO and FDA have different roles related to CIs and categorize CIs differently (ISO focuses on providing benchmark performance requirements for international standardization whereas the FDA focuses on the regulatory aspects of safety and efficacy for the U.S.). The FDA reviews all data, claims and labeling as part of its clearance process, while ISO compliance claims

are self-certified by the manufacturer. The FDA does not review ISO compliance claims. Understanding the roles and categorizations of ISO and FDA will help the SP team fully understand CI labeling and the information provided by each type of CI. ©

REFERENCES

1. International Organization for Standardization. *Sterilization of health care products — Vocabulary of terms used in sterilization and related equipment and process standards*. ISO 11139:2018. <https://www.iso.org/standard/66262.html>.
2. International Organization for Standardization. “Technical Committees.” <https://www.iso.org/technical-committees.html>.
3. International Organization for Standardization. *Sterilization of health care products — Chemical indicators—Part 1: General requirements*. ISO 11140-1:2014. <https://www.iso.org/standard/55080.html>.
4. U.S. Food and Drug Administration. *Guidance for Industry and FDA Staff: Premarket Notification [510(k)] Submissions for Chemical Indicators*. Document issued on December 19, 2003.



CRCST Self-Study Lesson Plan Quiz

Chemical Indicator Classifications: What Technicians Need to Know

Lesson No. CRCST 185 (Technical Continuing Education - TCE) · Lesson expires July 2025

1. Chemical indicators (CIs) provide information through a physical or chemical response to the sterilization process?
 - a. True
 - b. False
2. The ANSI/AAMI/ISO 11140-1 standard establishes:
 - a. Biological indicator (BI) classification
 - b. CI classification
 - c. The required use for all types of CIs and BIs
 - d. None of the above
3. CIs are intended to be used:
 - a. Inside packages
 - b. Outside packages
 - c. For special tests
 - d. All of the above
4. Which of the following organizations segment CIs into different classifications?
 - a. U.S. Food and Drug Administration (FDA) and the Association of periOperative Registered Nurses (AORN)
 - b. The Association for the Advancement of Medical Instrumentation (AAMI) and AORN
 - c. The FDA and the International Organization for Standardization (ISO)
 - d. AAMI and ISO
5. ISO 11140-1 has established how many types of CIs?
 - a. Three
 - b. Four
 - c. Five
 - d. Six
6. Type 6 indicators are used to:
 - a. Respond to one critical process
 - b. Respond to two critical processes
 - c. Respond to all critical processes
 - d. Test special functions
7. CIs can be used to:
 - a. Differentiate between processed and unprocessed items
 - b. Test special functions
 - c. Evaluate parameters and process variables at the point of placement
 - d. All of the above
8. A CI's stated value (SV) is:
 - a. The point where the CI shows the chamber exposure temperature
 - b. The manufacturer's statement of the CI's compatibilities
 - c. The point where the CI shows a pass result
 - d. The manufacturer's statement of the elements of a sterilization process that impact the lethality of the process
9. Which of the following directly affects the killing power of the sterilization process?
 - a. CI type
 - b. SV
 - c. Sterilization condition
 - d. Critical process variable
10. Type 1 indicators:
 - a. Have the same variables as Type 4 and 5 indicators
 - b. Require the use of SVs
 - c. Have a prescribed set of pass/fail conditions
 - d. Require a longer exposure time than Type 4 and 5 indicators
11. The FDA classifies which of the following?
 - a. Air removal tests
 - b. Process indicators
 - c. CIs
 - d. All of the above
12. Which FDA indicator type is used to distinguish between a processed and unprocessed package?
 - a. Air removal test
 - b. Process indicator
 - c. CI
 - d. BI
13. ISO and FDA have identical requirements for chemical indicators.
 - a. True
 - b. False
14. Requirements for CIs for VH2O2 sterilizers are covered under which section of the FDA guidance document?
 - a. Proper Use of VH2O2 Sterilizers
 - b. New Technology Sterilizers
 - c. Low-Temperature Sterilization
 - d. None of the above
15. Which organization reviews all claims and labeling of CIs before clearing the product?
 - a. ISO
 - b. ANSI/AAMI
 - c. FDA
 - d. All of the above

Request for paper/pencil scoring

(Please print or type information below)

- I have enclosed the scoring fee of \$15. (Please make checks payable to Purdue University. We regret that no refunds can be given.)
- Check here if you have a change of address.
- Check here if you wish to have your results emailed to you.

Detach quiz, fold, and return to:

Purdue Online – Central Service
Young Hall Room 405
155 S. Grant St.
West Lafayette, IN 47907

Name _____

Mailing Address (Be sure to include apartment number or post office box.) _____

City _____

State/Province _____

Zip Code/Postal Code _____

() Daytime Telephone _____

HSPA Membership Number _____

Email Address _____

If your name has changed in the last 12 months, please provide your former name. Purdue University is an equal access/equal opportunity institution.